Part Number	IFU-1392
Revision	05
Description	Instruction For Use, MSP, Multilanguage
Notes	N/A

Revision History

Rev	Effective Date	Summary of Change	DCO#	Author
00	5/6/14	Initial Release	14-047	M. Aldridge
01	6/12/15	Updated to reflect Warsaw address	15-018	M. Aldridge
02	12/28/15	Update revision to match IFU revision on label	15-110	S. Hartle
03	3/11/16	Added statement if product was re-used.	16-020	S. Hartle
04	8/24/17	Added new Emergo Address	17-082	S. Hartle
05	11/28/18	Added TL driver	18-110	C. Cook

Nextremity Solutions - Confidential Information

This document and the subject matter contained herein are the property of Nextremity Solutions and may not be copied, reproduced, or disclosed to others without the express written permission of a duly authorized agent of the company.

Sterile Implants & Instruments

Instructions for Use

MSP™ Metatarsal Shortening System

Please read this information carefully before using the MSP™ Metatarsal Shortening System

DESCRIPTION

The MSP™ Metatarsal Shortening System consists of a sterile multi component bone fixation implant (Metatarsal Shortening Plate, locking and cortical screws) and a set of accessory instruments designed for preparation of the implant site, and insertion of devices into the bone. These screws are all packaged individually and in separate packaging.

SYSTEM COMPONENTS

MSP™ Cortical Headed Screw:	Implantable Ti-6Al-4V ELI
MSP™ Locking Screw:	Implantable Ti-6Al-4V ELI
MSP™ Metatarsal Shortening Plate:	Implantable Ti-6Al-4V ELI
MSP™ Proximal Drill Guide:	Polycarbonate and surgical-grade stainless steel
MSP™ Drill Bushing:	Surgical-grade stainless steel
MSP™ Drill:	Surgical-grade stainless steel
MSP™ Driver:	Radel and surgical-grade stainless steel
MSP™ Driver, Torque Limiting:	Polyetherimide and Surgical Grade Stainless Steel

INDICATIONS FOR USE

The MSP™ Metatarsal Shortening System is indicated for fixation of fractures, osteotomies, nonunions, malunions and fusions of small bones and small bone segments such as the foot, particularly in osteopenic bone.

CONTRAINDICATIONS

1. Patient conditions including inadequate bone stock to receive implant.
2. Blood supply limitations and previous or active infections that may inhibit healing.
3. Surgical procedures other than for the indications listed.
4. Patients with conditions that limit their ability or willingness to follow postoperative care instructions.
5. Foreign-body sensitivity. Where material sensitivity is suspected, appropriate test should be made and sensitivity ruled out prior to implantation.

WARNINGS

This device should only be used by physicians with training in and a thorough understanding of orthopaedic surgery.

1. The MSP™ Metatarsal Shortening Plate is for SINGLE USE ONLY. DO NOT STERILIZE OR REUSE.
2. The MSP™ Metatarsal Shortening System has not been tested to withstand the forces needed for partial or full weight bearing or excessive activity until healing has occurred. The post-op regimen prescribed by the physician should be strictly followed to avoid stresses applied to the implant. Weight bearing post operatively should be at the discretion of the surgeon.
3. Do not use if the implants or instruments appear damaged upon inspection
4. Do not use in conjunction with implants of other metallic alloys
5. Any decision to remove the device should take into consideration the potential risk to the patient.
6. The MSP™ Metatarsal Shortening System has not been evaluated for safety and compatibility on the MR environment. The MSP™ Metatarsal Shortening System has not been tested for heating or migration in the MR environment.
7. The MSP™ Metatarsal Shortening System must be used only with provided insertion instruments.

CAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. A detailed surgical technique in print and/or electronic formats can be obtained by contacting a Nextremity Solutions sales representative.
2. Surgical instruments should only be used for their intended purpose.

POTENTIAL ADVERSE EFFECTS

1. Infection, deep and superficial
2. Foreign body reactions
3. Nonunion or delayed union which may lead to breakage of the implant
4. Bending or fracture of implant
5. Loosening or migration of the implant
6. Nerve damage due to surgical trauma
7. Inadequate healing
8. Pain, soft tissue discomfort or abnormal sensation due to the presence of the device

PREPARATION PRIOR TO USE

1. Inspect product for damage. Do not use if the product is damaged.
2. Users of this system are encouraged to contact their Nextremity Solutions representative if, in their professional judgment, they require a more comprehensive explanation of the surgical technique to be used with this device.
3. Insure that the patient is an appropriate candidate for the intended procedure based on the Indications, Contraindications, Warnings and Precautions listed above.
4. Weight bearing post operatively should be at the discretion of the surgeon.

DIRECTIONS FOR USE

WARNING: Use only the insertion instruments provided

1. Select the left or right implant kit based on the operative foot. The keel is oriented medially, and the dog-leg slot is oriented proximally. The distal edge of the plate is positioned approximately 5mm proximal to the metatarsal head. Thread the MSP™ Drill Bushing into the distal hole of the MSP™ Plate. Do not overtighten. Place the MSP™ Plate on the dorsal surface of the bone. Ensure the MSP™ Drill Bushing is in a vertical position. Apply dorsal and medial pressure to plate to ensure complete plate-to-bone contact while drilling. NOTE: Clamp may be used (not provided)
3. Use the 1.8mm drill bit to drill the distal hole. Drill a hole through the MSP™ Drill Bushing, into the bone. Use the markings on the Drill to measure the screw length at the top of the MSP™ Drill Bushing and select the appropriate length locking screw. Remove the drill bushing. Place the Locking Screw in the distal hole and fully tighten.
4. Thread the drill bushing into the second distal hole and use the drill bit to drill the second distal hole. Use the markings on the Drill to measure the screw length at the top of the MSP™ Drill Bushing and select the appropriate length locking screw. Remove drill bushing. Place a Locking Screw in the second distal hole and fully tighten.
5. NOTE: Ensure the keel is in contact with the medial cortex. Using the MSP™ Proximal Drill Guide, drill a hole into the bone in the proximal end of the dog-leg slot. Drill as far proximal as possible. Measure the screw length at the top of the MSP™ Proximal Drill Guide and select the appropriate length Cortical Screw. Place the Cortical Screw, but do not fully tighten. Note: Use the end with the single groove and apply medial pressure to the plate while drilling to ensure complete bone-to-plate contact. **ALTERNATIVE ORDER FOR SCREW PLACEMENT:** Drill and insert screws in the following order: 1. Most distal hole; 2. Dog-leg slot; 3. Second distal hole; 4. Proceed to step 6 as usual. NOTE: The alternate screw placement technique may help ensure consistent, optimal lateralization of the implant. For both the standard and alternate screw placement technique, it is critical to apply medial pressure to the plate while drilling to confirm complete bone-to-plate contact.
6. Using an oscillating saw, cut the bone obliquely through the cutting guide in the MSP™ Plate. The saw will be perpendicular to the dorsal surface of the MSP™ Plate.
7. Shorten the metatarsal by moving the bone proximally to the desired length. A total correction of 6mm is achieved if the screw reaches the most distal edge of the dog-leg slot.
8. Fully tighten the dog-leg screw.
9. Using the double-groove side of the MSP™ Proximal Drill Guide, insert the drill bit into the second proximal hole and drill. Use the markings on the Drill to measure the screw length at the top of the MSP™ Drill Bushing and select the appropriate length Cortical Screw. Place the Cortical Screw in the remaining proximal hole and fully tighten.

9.1 For supplemental compression of the osteotomy prior to drilling the final screw hole, a clamp may be used. The MSP™ Plate has a divot on the keel for this purpose. NOTE: Clamp not provided.

Optional device removal: In the event the device needs to be removed after placement has been completed.

1. Use the driver to carefully remove the screws from the bone.
2. Upon removal of the screws gently remove the MSP™ Metatarsal Shortening Plate from the bone with the use of small osteotomes.

CAUTION: IF THE DEVICE MUST BE REMOVED, DO NOT REUSE ANY PORTION OF THE DEVICE.

HOW SUPPLIED

The MSP™ Metatarsal Shortening System is available as a sterile surgical kit. This includes a plate and instruments, as well as a separate single sterile cortical headed and locking screws packaged individually. The components of the Metatarsal Shortening System are for single use only. DO NOT RE-STERILIZE. DO NOT REUSE. If product is re-used it may cause infection and may compromise the mechanical integrity of the product which may lead to failure.

The MSP™ Metatarsal Shortening System has not been evaluated for safety and compatibility in the MR environment. The MSP™ Metatarsal Shortening System has not been tested for heating or migration in the MR environment.

STORAGE AND HANDLING

Store in a cool, dry place and in a manner that protects the integrity of the packaging of all implants and instruments.

PACKAGING AND LABELING

1. Nextremity Solutions devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

Nextremity Solutions® and MSP™ are trademarks of Nextremity Solutions, Inc. These products are covered by one or more U.S. and international patents pending. All rights reserved. Printed in the USA

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